



A healthcare product



More efficiency in the  
microbiological control.

4% Sodium Hypochlorite



High-level disinfectant





## High-level disinfectant.

### Intended Use

#### Disinfection of:

- Semi-critical items (endoscopes, vaginal transducer, and others);
- Intermediate-level items (inhalation therapy material, ventilatory assistance, among others);
- Non-critical items (chamber pots, urinals, among others);
- Fixed surfaces (counters, stretchers, floors, walls, glasses, acrylic surfaces, among others);
- Touch screens, metal surfaces, among others;
- Critical areas (surgical theaters, ICU [intensive care units], hemodialysis, blood bank, among others);
- Semi-critical areas (recovery room, first aid station, among others);
- Non-critical areas (administrative area, general surfaces, restrooms, among others);
- Hospital cleaning areas;
- Nutrition and dietetics area;
- Specialty clinics (ophthalmology, gynecology, hemodialysis, etc.);
- Dental clinics;
- Ambulances;
- Clinical and pathological analysis labs;
- Forensic medicine institute;
- Others.

### How to use

While handling this product, please wear appropriate PPE, such as gloves, mask, goggles, apron, or gown.

- Add 5 ml of Atomic70 for each liter of water; Contacting time: 15 minutes;
- The solution can be used within a 24-hour period;
- Solution concentration is evidenced by measurement strips.

### Benefits

- Non-corrosive\*;
- Non-toxic\*;
- No need to demineralize the water;
- No need for rinsing it (except for semi-critical items);
- Eco-friendly;
- Eliminates biofilm;
- Best cost effectiveness;
- This product can be released directly in the sewage, not requiring neutralization.

\*At the recommended dose and contacting time

### Effective against many different Microorganisms:

- Staphylococcus aureus
- Salmonella choleraesuis
- Escherichia coli
- Pseudomonas aeruginosa
- Trichophyton mentagrophytes
- Candida albicans
- Mycobacterium smegmatis
- Mycobacterium massiliense
- Mycobacterium bovis
- Enterococcus faecium VRE (Vancomycin-resistant)
- Bacillus subtilis
- Clostridium sporogenes
- Acinetobacter baumannii
- klebsiella pneumoniae carbapenemase
- Enterococcus faecalis
- Staphylococcus Aureus MRSA
- Candida Auris
- Modified Vaccinia Virus Ankara
- Hepatitis B Virus
- Hepatitis C Virus
- Filoviridae (Ebola)
- Flavivirus (Chikungunya; Dengue)
- Human Immunodeficiency Virus
- Herpesviridae
- Human T - cell Leukemia Virus - HTLV - 1 Retrovirus
- Zika virus
- Influenza virus
- Coronavirus
- Paramyxoviridae (Measles)
- Rubella virus
- Burkholderia cepacia
- Listeria Monocytogenes



## Dioxide

Desenvolvendo Química Sustentável

ANVISA- M.S. No. 362940002  
Operation Permit by M.S.: 3.06.294-0

Technical Responsible Person:  
Alba Harumi Nishimura Saziki - CRQ/SP: 04405740

DIOXIDE INDÚSTRIA QUÍMICA LTDA.  
Avenida Vitória Rossi Martini, 641  
Distrito Industrial Vitória Martini  
Indaiatuba/SP - ZIP CODE: 13347-613  
Customer service: (19) 3936-1501 - [www.atomic70.com.br](http://www.atomic70.com.br)

In case of any emergencies, please call  
CEATOX (Toxicological Assistance Center): 0800 014 8110.