

## ATOMIC BRILHO CONCENTRATED ACID DETERGENT

# Revitalizes instruments

Removes stains and prevents darkening and oxidation



## HEALTHCARE PRODUCTS



Suitable for use on previously cleaned materials.



Designed to restore shine.



Removes mineral dirt and oxidation spots found on glassware, racks and surgical instruments made of stainless steel.



Used to neutralize residual alkalinity caused by the use of detergents in the reprocessing of these instruments.



**Dioxide**

Developing Sustainable Chemistry

[www.dioxide.com.br](http://www.dioxide.com.br)

  [dioxideindustriaquimica](https://www.instagram.com/dioxideindustriaquimica)



## ATOMIC BRILHO CONCENTRATED ACID DETERGENT




### When handling the product

Before starting use the product, assess the type of material on which the product will be applied in relation to:

- Reactivity;
- Resistance

- Chemical compatibility of each material.
- Instruments must be previously cleaned before being immersed in the solution.
- For more resistant stains, such as mineral plaque or oxidation, scrub the instruments with a soft brush (**do not use a hard-bristled metal brush**).
- Then rinse with plenty of water, completely removing any residue from the lumens and crevices of the instruments, and dry

### DIRECTIONS FOR USE

		
<b>Manual process at room temperature</b>	Dilute 50 mL of the product for each liter of water	Keep the instruments submerged for 20 minutes
<b>Manual and automated process at heated temperature (60oC to 80oC)</b>	Dilute 50 mL of product for every liter of water	Keep instruments submerged for 10 minutes
<b>After cleaning the instruments with alkaline detergent</b>	Dilute 4 mL of the product for each liter of water at room temperature or heated (60oC to 80oC)	Keep the solution circulating/ keep the instruments submerged for 2 to 5 minutes and rinse with plenty of water and dry



Improves joint performance



Extends shelf life of healthcare products (PPS)



High yield



Sanitizing product notified with ANVISA, No. 25351.329601/2023-47  
 Technician Responsible: Tatiane Cristina da Silva - CRQ/SP: 04476720

DIOXIDE INDÚSTRIA QUÍMICA LTDA  
 Avenida Vitória Rossi Martini, 641 - CEP:13347-613  
 Distrito Industrial Vitória Martini - Indaiatuba/SP  
 SAC (19) 3936-1501 - www.dioxide.com.br  
 In case of emergency, call CEATOX: 0800-722-8001

