

#### **CLEANS AND DISINFECTS HEALTHCARE ENVIRONMENTS**



## WHAT ARE THE ADVANTAGES OF ATOMIC READY?







**FAST ACTING** 



BIODEGRADABLE



**ECONOMIC** 









### **Characteristics**

- ✓ NEUTRAL pH
- ✓ ALCOHOL-FREE
- ✓ POLYMERIC BIGUANIDE
- ✓ NON-IONIC SURFACTANT
- √ 5<sup>th</sup> GENERATION AMMONIA QUATERNARY

- ✓ NO RINSING
- **✓** NO DYES
- **✓** NO PERFUME
- **✓** NON-CORROSIVE
- **✓** NON-FLAMMABLE

# Atomic Ready is recommended for:



Hospitals



Operating room



Specialized clinics



Noncritical articles



<u>Ambulances</u>



Outpatient department



#### How to use:

To clean and disinfect fixed surfaces, apply Atomic Ready directly on the surface with a spray, pulverizer or a clean cloth wet made moist this product. To clean and disinfect noncritical articles, wash using Atomic Ready and let the product act for 2 minutes, then rinse and dry the articles. \*Wear gloves during usage.

PRODUCED AND DISTRIBUTED BY THE BRAZILIAN INDUSTRY: DIOXIDE INDÚSTRIA QUÍMICA LTDA. BRAZILIAN IND. / ANVISA | H.M. REGISTER No. 362940005 | CSH (19) 3936-1501





